

For Research Use Only

MKNK1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86966-3-PBS



Basic Information

Catalog Number:

86966-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0184

GenBank Accession Number:

BC002755

GeneID (NCBI):

8569

UNIPROT ID:

Q9BUB5

Full Name:

MAP kinase interacting
serine/threonine kinase 1

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

252096C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MNK1 concentration was determined to be 8.4 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 9.4 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 2.5 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02797-1, MNK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86966-3-PBS. Detection antibody: 86966-2-PBS. Standard: Ag0184. Range: 0.312-40 ng/mL.