

For Research Use Only

TIP30 Recombinant monoclonal antibody

Catalog Number: 86972-1-RR



Basic Information

Catalog Number: 86972-1-RR	GenBank Accession Number: BC002439	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10553	CloneNo.: 252091A12
Isotype: IgG	UNIPROT ID: Q9BUP3	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: AG0353	Full Name: HIV-1 Tat interactive protein 2, 30kDa	
	Calculated MW: 27 kDa	
	Observed MW: 28 kDa, 32 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, A431 cells, A549 cells, SW 1990 cells, A-673 cells
Species Specificity: human	

Background Information

TIP30, namely HTATIP2 or CC3, was firstly discovered in small cell lung carcinoma, using differential analyses between highly metastatic human variant cells and less metastatic classic cells. TIP30 was found to be downregulated in various tumors and is considered to be a tumor suppressor due to its pro-apoptotic activity and its anti-metastatic and anti-angiogenic capacities.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86972-1-RR (TIP30 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.