

For Research Use Only

DNAJC5B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86975-5-PBS



Basic Information

Catalog Number:

86975-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10744

GenBank Accession Number:

BC016742

GeneID (NCBI):

85479

UNIPROT ID:

Q9UF47

Full Name:

DnaJ (Hsp40) homolog, subfamily C, member 5 beta

Calculated MW:

199 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

252079F12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DNAJC5B concentration was determined to be 137.7 ng/mL in mouse testis tissue extract based on a 4.2 mg/mL extract load, mouse lung tissue extract were served as a negative control.



Cytometric bead array standard curve of MP02805-1, DNAJC5B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86975-2-PBS. Detection antibody: 86975-5-PBS. Standard: Ag10744. Range: 0.312-40 ng/mL.