

For Research Use Only

Villin Recombinant monoclonal antibody

Catalog Number: 86978-3-RR



Basic Information

Catalog Number: 86978-3-RR	GenBank Accession Number: BC017303	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 7429	CloneNo.: 251823E12
Isotype: IgG	UNIPROT ID: P09327	Recommended Dilutions: WB: 1:5000-1:50000
Immunogen Catalog Number: AG9610	Full Name: villin 1	
	Calculated MW: 827aa, 93 kDa; 826aa, 93 kDa	
	Observed MW: 93-95 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: human, mouse, rat	WB : mouse colon tissue, rat colon tissue, mouse kidney tissue, rat kidney tissue

Background Information

Villin 1 (VIL1) is a 95-kd F-actin bundling and severing protein and its expression is restricted to epithelial cells with a brush border, like epithelial cells of the intestinal mucosa, gall bladder, renal proximal tubules and ductuli efferentes of the testis. VIL1 has been reported to be an epithelial cell-specific anti-apoptotic protein, and to have an important function in regulating actin dynamics, cell morphology, epithelial-to-mesenchymal transitions, cell migration and cell survival. In addition, VIL1 is a useful diagnostic marker for various cancer, like cervical and endometrial adenocarcinomas, renal cell carcinoma. VIL1 was recently identified as a novel biomarker predictive for postoperative recurrence and poorer prognosis of high serum AFP associated HCC.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86978-3-RR (VIL1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.