

For Research Use Only

Rat CD43 Recombinant monoclonal antibody, PBS Only

Catalog Number:86985-1-PBS



Basic Information

Catalog Number: 86985-1-PBS	GenBank Accession Number: P13838	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 24796	CloneNo.: 251932G11
Isotype: IgG	UNIPROT ID: P13838	
Immunogen Catalog Number: EG5149	Full Name: sialophorin	
	Calculated MW: 38 kDa	
	Observed MW: 90-110 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
rat

Background Information

CD43, also known as leukosialin, sialophorin, leukocyte sialoglycoprotein, SPN and W3/13 antigen, is a mucin-like type I transmembrane protein. Rat CD43 is a sialoglycoprotein with an apparent molecular weight of ~100 kDa and is found on thymocytes, T lymphocytes, neutrophils, plasma cells and myelomas but not on B lymphocytes or most other tissues (PMID: 2965006). CD43 is involved in the contact between cells participating in a series of events such as signaling, adherence, and host-parasite interactions (PMID: 10580634; 31197200).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat thymus tissue was subjected to SDS PAGE followed by western blot with 86985-1-RR (CD43 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86985-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86985-1-RR against Rat CD43 were performed. The affinity constant is 0.256 nM.