

For Research Use Only

CA5B Recombinant monoclonal antibody, PBS Only

Catalog Number:86996-1-PBS



Basic Information

Catalog Number: 86996-1-PBS	GenBank Accession Number: BC028142	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 11238	CloneNo.: 251915F10
Isotype: IgG	UNIPROT ID: Q9Y2D0	
Immunogen Catalog Number: AG4166	Full Name: carbonic anhydrase VB, mitochondrial	
	Calculated MW: 317 aa, 36 kDa	
	Observed MW: 32-35 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

CA5B(Carbonic anhydrase 5B, mitochondrial) is also named as CA-VB and belongs to the alpha-carbonic anhydrase family.CA5B encodes a 317-amino acid protein that has 64% sequence similarity with CA5 (114761), a mitochondrial CA, and contains a 33-amino-acid residue signal sequence. Western blot analysis revealed a 36-kD protein precursor and a 32-kD mature protein (PMID:10409679).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86996-1-RR (CA5B antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86996-1-PBS in a different storage buffer formulation.