

For Research Use Only

PPARGC1B Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86997-3-PBS



Basic Information

Catalog Number:

86997-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17566

GenBank Accession Number:

BC132971

GeneID (NCBI):

133522

UNIPROT ID:

Q86YN6

Full Name:

peroxisome proliferator-activated receptor gamma, coactivator 1 beta

Calculated MW:

1023 aa, 113 kDa

Purification Method:

Protein A purification

CloneNo.:

252087B12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PPARGC1B concentration was determined to be 0.85 ng/mL in HeLa cell extract based on a 1.20 mg/mL extract load, 1.02 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load and 0.45 ng/mL in SH-SY5Y cell extract based on a 1.40 mg/mL extract load.



Sandwich ELISA standard curve of MP02897-1, Human PPARGC1B Recombinant Matched Antibody Pair - PBS only. 86997-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17566. 86997-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL