

For Research Use Only

CCL23 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86998-2-PBS



Basic Information

Catalog Number: 86998-2-PBS	GenBank Accession Number: BC093970	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 6368	CloneNo.: 250795A1
Isotype: IgG	UNIPROT ID: P55773	
Immunogen Catalog Number: EG4244	Full Name: chemokine (C-C motif) ligand 23	
	Calculated MW: 120 aa, 13 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CCL23 concentration was determined to be 537.2 pg/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02815-1, Human CCL23 Recombinant Matched Antibody Pair - PBS only. 86998-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4244. 86998-1-PBS was HRP conjugated as the detection antibody. Range: 15.6-1000 pg/mL.