

For Research Use Only

PDCD5 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 86999-2-PBS



Basic Information

Catalog Number:

86999-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3122

GenBank Accession Number:

BC015519

GeneID (NCBI):

9141

UNIPROT ID:

O14737

Full Name:

programmed cell death 5

Calculated MW:

125 aa, 14 kDa

Purification Method:

Protein A purification

CloneNo.:

251624B5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of twelve individual healthy human donors was measured. The PDCD5 concentration of detected samples was determined to be 133.4 pg/mL with a range of 73.5-195.5 pg/mL.



Sandwich ELISA standard curve of MP02813-1, Human PDCD5 Recombinant Matched Antibody Pair - PBS only. 86999-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3122. 86999-1-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.