

For Research Use Only

BEND3 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87001-3-PBS



Basic Information

Catalog Number: 87001-3-PBS	GenBank Accession Number: BC131628	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 57673	CloneNo.: 252025F2
Isotype: IgG	UNIPROT ID: Q5T5X7	
Immunogen Catalog Number: AG19137	Full Name: BEN domain containing 3	
	Calculated MW: 828 aa, 94 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BEND3 concentration was determined to be 76.4 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 18.2 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP03045-1, BEND3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87001-1-PBS. Detection antibody: 87001-3-PBS. Standard: Ag19137. Range: 1.562-200 ng/mL