

For Research Use Only

PRKCSH Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:87004-2-PBS



Basic Information

Catalog Number: 87004-2-PBS	GenBank Accession Number: BC013586	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5589	CloneNo.: 252090A9
Isotype: IgG	UNIPROT ID: P14314	
Immunogen Catalog Number: AG2796	Full Name: protein kinase C substrate 80K-H	
	Calculated MW: 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PRKCSH concentration was determined to be 34.6 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load, 31.1 ng/mL in mouse liver tissue extract based on a 8.4 mg/mL extract load, 6.5 ng/mL in HCT 116 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02821-1, PRKCSH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87004-2-PBS. Detection antibody: 87004-1-PBS. Standard: Ag2796. Range: 0.469-60 ng/mL.