

For Research Use Only

Vitronectin Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87005-2-PBS



Basic Information

Catalog Number:

87005-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4682

GenBank Accession Number:

NM_000638.4

GeneID (NCBI):

7448

UNIPROT ID:

P04004

Full Name:

vitronectin

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

252166C10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Vitronectin/S-protein concentration was determined to be 430.3 ug/mL in human plasma.



Cytometric bead array standard curve of MP02822-1, Vitronectin/S-protein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87005-2-PBS. Detection antibody: 87005-1-PBS. Standard: Eg4682. Range: 0.625-80 ng/mL.