

For Research Use Only

ADAM8 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87008-2-PBS



Basic Information

Catalog Number:

87008-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2912

GenBank Accession Number:

NM_001109.4

GeneID (NCBI):

101

UNIPROT ID:

P78325-1

Full Name:

ADAM metallopeptidase domain 8

Calculated MW:

89 kDa

Purification Method:

Protein A purification

CloneNo.:

251352A5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Serum of eight individual healthy human donors was measured. The human ADAM8 concentration of detected samples was determined to be 1.88 ng/mL with a range of 0.36 - 3.20 ng/mL.



Sandwich ELISA standard curve of MP02833-1, Human ADAM8 Recombinant Matched Antibody Pair - PBS only. 87008-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2912. 87008-2-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL