

For Research Use Only

ELK1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87009-3-PBS



Basic Information

Catalog Number: 87009-3-PBS	GenBank Accession Number: BC056150	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2002	CloneNo.: 252085A8
Isotype: IgG	UNIPROT ID: P19419	
Immunogen Catalog Number: AG26683	Full Name: ELK1, member of ETS oncogene family	
	Calculated MW: 45 kDa	
	Observed MW: 58-60 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ELK1 is a nuclear transcription factor of the ETS-domain family and the ternary complex factor (TCF) sub-family. It binds to purine-rich DNA sequences and forms a ternary complex with serum response factor (SRF) on the serum response element (SRE) located in the promoter regions of immediate early genes such as FOS and IER2. Upon activation of the JNK or MAPK/ERK signaling cascades, ELK1 is phosphorylated and rapidly induces target-gene transcription, thereby regulating cell proliferation, differentiation, and survival.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87009-3-RR (ELK1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87009-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87009-3-RR against Human ELK1 were performed. The affinity constant is 4.94 nM.