

For Research Use Only

# PSMB9 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87013-3-PBS



## Basic Information

Catalog Number:

87013-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6048

GenBank Accession Number:

BC065513

GeneID (NCBI):

5698

UNIPROT ID:

P28065

Full Name:

proteasome (prosome, macropain)  
subunit, beta type, 9 (large  
multifunctional peptidase 2)

Calculated MW:

23 kDa

Purification Method:

Protein A purification

CloneNo.:

252148A9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PSMB9 concentration was determined to be 5.1 ng/mL in HL-60 cell extract based on a 1.3 mg/mL extract load, HEK-293T cell extract were served as a negative control.



Cytometric bead array standard curve of MP02846-1, PSMB9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87013-1-PBS. Detection antibody: 87013-3-PBS. Standard: Ag6048. Range: 0.312-40 ng/mL.