

For Research Use Only

P130 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87016-3-PBS



Basic Information

Catalog Number:

87016-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2382

GenBank Accession Number:

BC006769

GeneID (NCBI):

9221

UNIPROT ID:

Q14978

Full Name:

nucleolar and coiled-body phosphoprotein 1

Calculated MW:

73 kDa

Purification Method:

Protein A purification

CloneNo.:

252136A10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean P130 concentration was determined to be 266.96 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load and 164.09 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02870-1, Human P130 Recombinant Matched Antibody Pair - PBS only. 87016-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2382. 87016-2-PBS was HRP conjugated as the detection antibody. Range: 0.39-25 ng/mL