

For Research Use Only

ZCCHC8 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87017-4-PBS



Basic Information

Catalog Number:

87017-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19985

GenBank Accession Number:

BC065918

GeneID (NCBI):

55596

UNIPROT ID:

Q6NZY4

Full Name:

zinc finger, CCHC domain containing
8

Calculated MW:

707 aa, 79 kDa

Purification Method:

Protein A purification

CloneNo.:

252120G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ZCCHC8 concentration was determined to be 46.0 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 22.1 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 16.3 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02824-1, ZCCHC8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87017-4-PBS. Detection antibody: 87017-3-PBS. Standard: Ag19985. Range: 0.156-20 ng/mL.