

For Research Use Only

BMX Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87023-2-PBS



Basic Information

Catalog Number:

87023-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26407

GenBank Accession Number:

BC016652

GeneID (NCBI):

660

UNIPROT ID:

P51813

Full Name:

BMX non-receptor tyrosine kinase

Calculated MW:

675 aa, 78 kDa

Purification Method:

Protein A purification

CloneNo.:

252193D1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BMX concentration was determined to be 35.71 ng/mL in HT-29 cell extract based on a 1.5 mg/mL extract load.



Sandwich ELISA standard curve of MP02871-1, Human BMX Recombinant Matched Antibody Pair - PBS only. 87023-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag26407. 87023-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL