

For Research Use Only

BRD3 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87030-4-PBS



Basic Information

Catalog Number:

87030-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2433

GenBank Accession Number:

BC032124

GeneID (NCBI):

8019

UNIPROT ID:

Q15059

Full Name:

bromodomain containing 3

Calculated MW:

556 aa, 61 kDa

Purification Method:

Protein A purification

CloneNo.:

252143B6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean BRD3 concentration was determined to be 16.21 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load and 8.56 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02845-2, Human BRD3 Recombinant Matched Antibody Pair - PBS only. 87030-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2433. 87030-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL