

For Research Use Only

NUSAP1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87034-1-PBS



Basic Information

Catalog Number:

87034-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2654

GenBank Accession Number:

BC001308

GeneID (NCBI):

51203

UNIPROT ID:

Q9BX56

Full Name:

nucleolar and spindle associated protein 1

Calculated MW:

440 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

252131D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NUSAP1 concentration was determined to be 99.0 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 84.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 38.4 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02838-1, NUSAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87034-2-PBS. Detection antibody: 87034-1-PBS. Standard: Ag2654. Range: 0.625-80 ng/mL.