

For Research Use Only

NUSAP1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87034-5-PBS



Basic Information

Catalog Number:

87034-5-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2654

GenBank Accession Number:

BC001308

GeneID (NCBI):

51203

UNIPROT ID:

Q9BXS6

Full Name:

nucleolar and spindle associated protein 1

Calculated MW:

440 aa, 49 kDa

Purification Method:

Protein A purification

CloneNo.:

252131D7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean NUSAP1 concentration was determined to be 9.4 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 19.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load and 0.2 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02838-2, Human NUSAP1 Recombinant Matched Antibody Pair - PBS only. 87034-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2654. 87034-4-PBS was HRP conjugated as the detection antibody. Range: 46.9-1500 pg/mL.