

For Research Use Only

TRPM2 Recombinant monoclonal antibody

Catalog Number:87036-1-RR



Basic Information

Catalog Number:

87036-1-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC041570

GeneID (NCBI):

7226

UNIPROT ID:

O94759

Full Name:

transient receptor potential cation channel, subfamily M, member 2

Calculated MW:

171 kDa

Observed MW:

171 kDa

Purification Method:

Protein A purification

CloneNo.:

251122F2

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MDA-MB-231 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MDA-MB-231 cells were subjected to SDS PAGE followed by western blot with 87036-1-RR (TRPM2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.