

For Research Use Only

MYD88 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87040-3-PBS



Basic Information

Catalog Number:

87040-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19770

GenBank Accession Number:

BC013589

GeneID (NCBI):

4615

UNIPROT ID:

Q99836

Full Name:

myeloid differentiation primary response gene (88)

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

252132E3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MYD88 concentration was determined to be 907.15 ng/mL in Jurkat cell extract based on a 1.20 mg/mL extract load, 2,599.63 ng/mL in Raji cell extract based on a 1.20 mg/mL extract load and 123.62 ng/mL in HCT 116 cell extract based on a 1.30 mg/mL extract load.



Sandwich ELISA standard curve of MP03184-1, Human MYD88 Recombinant Matched Antibody Pair - PBS only. 87040-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19770. 87040-3-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL