

For Research Use Only

DNAJA2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87041-1-PBS



Basic Information

Catalog Number: 87041-1-PBS	GenBank Accession Number: BC015809	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10294	CloneNo.: 252019B6
Isotype: IgG	UNIPROT ID: O60884	
Immunogen Catalog Number: AG2877	Full Name: DnaJ (Hsp40) homolog, subfamily A, member 2	
	Calculated MW: 417 aa, 46 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DNAJA2 concentration was determined to be 317.8 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load, 201.1 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 118.7 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02847-1, DNAJA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87041-2-PBS. Detection antibody: 87041-1-PBS. Standard: Ag2877. Range: 1.562-200 ng/mL.