

For Research Use Only

# DNAJA2 Recombinant monoclonal antibody

Catalog Number:87041-2-RR



## Basic Information

Catalog Number:

87041-2-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2877

GenBank Accession Number:

BC015809

GeneID (NCBI):

10294

UNIPROT ID:

O60884

Full Name:

DnaJ (Hsp40) homolog, subfamily A, member 2

Calculated MW:

417 aa, 46 kDa

Observed MW:

46 kDa

Purification Method:

Protein A purification

CloneNo.:

252019B10

Recommended Dilutions:

WB: 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB: A375 cells, MCF-7 cells, HepG2 cells, Jurkat cells, K-562 cells, fetal human brain tissue, mouse kidney tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87041-2-RR (DNAJA2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.