

For Research Use Only

# SAP30BP Recombinant monoclonal antibody

Catalog Number:87054-1-RR



## Basic Information

<b>Catalog Number:</b> 87054-1-RR	<b>GenBank Accession Number:</b> BC007592	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 29115	<b>CloneNo.:</b> 252137B1
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> Q9UHR5	<b>Recommended Dilutions:</b> WB: 1:5000-1:50000
<b>Immunogen Catalog Number:</b> AG15051	<b>Full Name:</b> SAP30 binding protein	
	<b>Calculated MW:</b> 308 aa, 34 kDa	
	<b>Observed MW:</b> 45 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : MCF-7 cells, A549 cells, Jurkat cells, HEK-293T cells, mouse spleen tissue, rat spleen tissue, mouse kidney tissue
<b>Species Specificity:</b> human, mouse, rat	

## Background Information

SAP30BP, also named as SAP30-binding protein, is a 308 amino acid protein, which belongs to the HCNGP family. SAP30BP localizes in Nucleus and Interacts with SAP30. DDX54 Induces cell death and may act as a transcriptional corepressor of a gene related to cell survival. SAP30BP may be involved in the regulation of beta-2-microglobulin genes.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
*Aliquoting is unnecessary for -20°C storage*

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87054-1-RR (SAP30BP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 87054-1-RR against Human SAP30BP were performed. The affinity constant is 1.23 nM.