

For Research Use Only



Integrin alpha-1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87056-1-PBS

Basic Information

Catalog Number: 87056-1-PBS	GenBank Accession Number: BC137121	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3672	CloneNo.: 251984E2
Isotype: IgG	UNIPROT ID: P56199	
Immunogen Catalog Number: EG2089	Full Name: integrin, alpha 1	
	Calculated MW: 1179 aa, 131 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Integrin alpha-1 concentration was determined to be 51.4 ng/mL in human plasma.



Cytometric bead array standard curve of MP02854-1, Integrin alpha-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87056-2-PBS. Detection antibody: 87056-1-PBS. Standard: Eg2089. Range: 0.781-100 ng/mL.