

For Research Use Only

UBL3 Recombinant monoclonal antibody, PBS Only

Catalog Number:87070-3-PBS



Basic Information

Catalog Number: 87070-3-PBS	GenBank Accession Number: NM_007106.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5412	CloneNo.: 252163G4
Isotype: IgG	UNIPROT ID: O95164	
Immunogen Catalog Number: EG5932	Full Name: ubiquitin-like 3	
	Calculated MW: 13 kDa	
	Observed MW: 13 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Ubiquitin-like 3 (UBL3)/membrane-anchored Ub-fold protein (MUB) acts as a post-translational modification (PTM) factor to regulate efficient protein sorting to sEVs (small extracellular vesicles). (PMID: 31363817, PMID: 30258067).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87070-3-RR (UBL3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87070-3-PBS in a different storage buffer formulation.