

For Research Use Only

CD21 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87097-2-PBS



Basic Information

Catalog Number:

87097-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2980

GenBank Accession Number:

NM_001877.5

GeneID (NCBI):

1380

ENSEMBL Gene ID:

ENSG00000117322

UNIPROT ID:

P20023-1

Full Name:

complement component (3d/Epstein Barr virus) receptor 2

Calculated MW:

113 kDa

Purification Method:

Protein A purification

CloneNo.:

252181B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD21 concentration was determined to be 131.8 ng/mL in human plasma.



Cytometric bead array standard curve of MP02864-1, CD21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87097-3-PBS. Detection antibody: 87097-2-PBS. Standard: Eg2980. Range: 0.938-60 ng/mL.