

For Research Use Only

# Apolipoprotein CIII Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87100-3-PBS



## Basic Information

<b>Catalog Number:</b> 87100-3-PBS	<b>GenBank Accession Number:</b> NM_000040	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 345	<b>CloneNo.:</b> 252153D7
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P02656	
<b>Immunogen Catalog Number:</b> EG3387	<b>Full Name:</b> apolipoprotein C-III	
	<b>Calculated MW:</b> 11 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean Apolipoprotein CIII concentration was determined to be 218.5 mg/mL in human serum.



Cytometric bead array standard curve of MP02866-1, Apolipoprotein CIII Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87100-3-PBS. Detection antibody: 87100-2-PBS. Standard: Eg3387. Range: 0.625-40 ng/mL.