

For Research Use Only

MUC13 Recombinant monoclonal antibody

Catalog Number: 87107-2-RR



Basic Information

Catalog Number: 87107-2-RR	GenBank Accession Number: NM_033049.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 56667	CloneNo.: 252211G3
Isotype: IgG	UNIPROT ID: Q9H3R2	Recommended Dilutions: WB: 1:1000-1:4000 IHC: 1:500-1:2000 IF/ICC: 1:250-1:1000
Immunogen Catalog Number: EG6586	Full Name: mucin 13, cell surface associated	
	Calculated MW: 55 kDa	
	Observed MW: 120 kDa	

Applications

Tested Applications:
WB, IHC, IF/ICC, ELISA

Species Specificity:
human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HT-29 cells,
IHC : human colon tissue,
IF/ICC : HT-29 cells,

Background Information

Mucins (MUC) are highly glycosylated proteins that are abundantly expressed in mammalian epithelial cells. Mucin 13 (MUC13) is a member of the transmembrane glycoprotein mucins, which are normally expressed in the large intestine, trachea, kidney, small intestine, and gastric epithelium (PMID: 22027689). MUC13 is involved in the tumorigenesis of multiple malignancies, including pancreatic cancer, colorectal cancer, renal cancer, and lung cancer (PMID: 22027689; PMID: 31427737; PMID: 28205224; PMID: 34966453).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HT-29 cells were subjected to SDS PAGE followed by western blot with 87107-2-RR (MUC13 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 87107-2-RR (MUC13 antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 87107-2-RR (MUC13 antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed HT-29 cells using MUC13 antibody (87107-2-RR, Clone: 25221G3) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



Biolayer interferometry (BLI) kinetic assays of 87107-1-RR against Human MUC13 were performed. The affinity constant is below 1 pM.