

For Research Use Only

ART3 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87108-1-PBS



Basic Information

Catalog Number: 87108-1-PBS	GenBank Accession Number: BC008397	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 419	CloneNo.: 252232A3
Isotype: IgG	UNIPROT ID: Q13508	
Immunogen Catalog Number: AG8775	Full Name: ADP-ribosyltransferase 3	
	Calculated MW: 389 aa, 44 kDa	
	Observed MW: 38-47 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Background Information

ART3 is a GPI-anchored mono-ADP-ribosyltransferase enriched on the sperm cell surface, where it participates in spermatid maturation and membrane remodeling. Localized to lipid rafts, ART3 likely modulates cell-cell signaling and extracellular protein modification. Its upregulation in multiple cancers suggests additional roles in migration, invasion, and tumor aggressiveness. Despite incomplete characterization of its enzymatic activity, ART3 remains a key regulator at the interface of reproductive biology, membrane signaling, and pathology.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 87108-1-RR (ART3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87108-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87108-1-RR against Human ART3 were performed. The affinity constant is 0.205 nM.