

For Research Use Only

APLP1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87113-1-PBS



Basic Information

Catalog Number: 87113-1-PBS	GenBank Accession Number: NM_005166.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 333	CloneNo.: 252210C4
Isotype: IgG	UNIPROT ID: P51693-1	
Immunogen Catalog Number: EG5634	Full Name: amyloid beta (A4) precursor-like protein 1	
	Calculated MW: 72 kDa	
	Observed MW: 80 kDa	

Applications

Tested Applications:
WB, IHC, IF-P, Indirect ELISA

Species Specificity:
human, mouse

Background Information

APLP1 (Amyloid Precursor-Like Protein 1) is a member of the amyloid precursor protein (APP) family, which also includes APP and APLP2. It is a type I transmembrane glycoprotein that is highly and predominantly expressed in the central nervous system, especially at the postsynaptic density of neurons. Functionally, APLP1 plays important roles in synaptic adhesion, synaptogenesis, and neurite outgrowth. It can form homophilic trans-interactions between neighboring cells, contributing to cell-cell adhesion and the formation of synaptic structures. Its proteolytic processing by α -, β -, and γ -secretases releases intracellular fragments that may participate in transcriptional regulation and signal transduction, potentially influencing gene expression and neuronal survival.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




fetal human brain tissue were subjected to SDS PAGE followed by western blot with 87113-1-RR (APLP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87113-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 87113-1-RR (APLP1 antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 87113-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 87113-1-RR (APLP1 antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 87113-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse brain tissue using APLP1 antibody (87113-1-RR, Clone: 252210C4) at dilution of 1:200 and Coralite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 87113-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87113-1-RR against Human APLP1 were performed. The affinity constant is below 1 pM.