

For Research Use Only

RBL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87119-1-PBS



Basic Information

Catalog Number: 87119-1-PBS	GenBank Accession Number: BC032247	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5933	CloneNo.: 252269A11
Isotype: IgG	UNIPROT ID: P28749	
Immunogen Catalog Number: AG4146	Full Name: retinoblastoma-like 1 (p107)	
	Calculated MW: 1068 aa, 119 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean p107 concentration was determined to be 50.3 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 29.2 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02868-1, p107 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87119-2-PBS. Detection antibody: 87119-1-PBS. Standard: Ag4146. Range: 0.469-60 ng/mL.