

For Research Use Only

Mouse Flt-4 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87120-2-PBS



Basic Information

| | | |
|--|---|---|
| Catalog Number: 87120-2-PBS | GenBank Accession Number: NM_008029.3 | Purification Method: Protein A purification |
| Source: Rabbit | GeneID (NCBI): 14257 | CloneNo.: 252201B7 |
| Isotype: IgG | UNIPROT ID: P35917 | |
| Immunogen Catalog Number: EG6355 | Full Name: FMS-like tyrosine kinase 4 | |
| | Calculated MW: 153 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FLT4/VEGFR3 concentration was determined to be 741.1 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02875-1, MOUSE FLT4/VEGFR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87120-2-PBS. Detection antibody: 87120-1-PBS. Standard: Eg6355. Range: 0.469-60 ng/mL.