

For Research Use Only

# PSMC3 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87128-2-PBS



## Basic Information

<b>Catalog Number:</b> 87128-2-PBS	<b>GenBank Accession Number:</b> BC008713	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 5702	<b>CloneNo.:</b> 252375C2
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P17980	
<b>Immunogen Catalog Number:</b> AG21408	<b>Full Name:</b> proteasome (prosome, macropain) 26S subunit, ATPase, 3	
	<b>Calculated MW:</b> 49 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PSMC3 concentration was determined to be 81.4 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 33.8 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 28.3 ng/mL in HCT 116 cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02881-1, PSMC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87128-2-PBS. Detection antibody: 87128-1-PBS. Standard: Ag21408. Range: 0.938-60 ng/mL.