

For Research Use Only

TSPAN14 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87131-3-PBS



Basic Information

Catalog Number: 87131-3-PBS	GenBank Accession Number: BC002920	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 81619	CloneNo.: 251354E9
Isotype: IgG	UNIPROT ID: Q8NG11	
	Full Name: tetraspanin 14	
	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TSPAN14 concentration was determined to be 217.2 pg/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02883-1, TSPAN14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87131-4-PBS. Detection antibody: 87131-3-PBS. Standard: SY02596. Range: 39-2500 pg/mL