

For Research Use Only

TSPAN14 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87131-4-PBS



Basic Information

Catalog Number:

87131-4-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC002920

GeneID (NCBI):

81619

UNIPROT ID:

Q8NG11

Full Name:

tetraspanin 14

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

251354A10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TSPAN14 concentration was determined to be 217.2 pg/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02883-1, TSPAN14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87131-4-PBS. Detection antibody: 87131-3-PBS. Standard: SY02596. Range: 39-2500 pg/mL