

For Research Use Only

IFIH1/MDA5 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87135-1-PBS



Basic Information

Catalog Number: 87135-1-PBS	GenBank Accession Number: BC111750	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 64135	CloneNo.: 252265E9
Isotype: IgG	UNIPROT ID: Q9BYX4	
Immunogen Catalog Number: AG16467	Full Name: interferon induced with helicase C domain 1	
	Calculated MW: 1025 aa, 117 kDa	
	Observed MW: 135 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

IFIH1 (Interferon-induced helicase C domain-containing protein 1) is a putative RNA helicase that is upregulated in response to treatment with IFNB or IFNB and MEZ. It is also named as MDA5 and RH116. Ectopic expression of MDA5 in melanoma cells resulted in reduced colony formation, suggesting an interaction of the CARD and apoptotic signal molecules. Functional analysis indicated that MDA5 is an RNA-dependent ATPase. IFIH1 has an apparent molecular mass of 117-130 kDa, and always other bands (70 kDa and 90 kDa) can be detected as cleaved products (PMID: 17267501). This antibody is specific to IFIH1.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87135-1-RR (IFIH1/MDA5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87135-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87135-1-RR against Human IFIH1 were performed. The affinity constant is 5.16 nM.