

For Research Use Only

SPTLC1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87141-1-PBS



Basic Information

Catalog Number: 87141-1-PBS	GenBank Accession Number: BC007085	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 10558	CloneNo.: 252268G2
Isotype: IgG	UNIPROT ID: O15269	
Immunogen Catalog Number: AG1162	Full Name: serine palmitoyltransferase, long chain base subunit 1	
	Calculated MW: 53 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

SPTLC1 is a subunit of serine palmitoyltransferase (SPT) which is the key enzyme in sphingolipid biosynthesis and is essential for embryogenesis and cell survival. Mutations in the SPTLC1 gene (C133W, C133Y, V144D, and G387A) were reported to be responsible for the development of an inherited sensory neuropathy (hereditary sensory neuropathy type I, HSN1) (PMID: 39959268). Pathogenic variants in SPTLC1 are causative for hereditary sensory and autonomic neuropathy, juvenile amyotrophic lateral sclerosis, macular telangiectasia type 2, or Charcot-Marie-Tooth disease (PMID: 31751474). Western blot analysis detected a specific band at ~53 kDa (PMID: 36197001).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87141-1-RR (SPTLC1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87141-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87141-1-RR against Human SPTLC1 were performed. The affinity constant is 0.487 nM.