

For Research Use Only

MUM1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 87147-1-PBS



Basic Information

Catalog Number: 87147-1-PBS	GenBank Accession Number: BC008098	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 84939	CloneNo.: 252371H1
Isotype: IgG	UNIPROT ID: Q2TAK8	
Immunogen Catalog Number: AG3371	Full Name: melanoma associated antigen (mutated) 1	
	Calculated MW: 18 kDa, 79 kDa	
	Observed MW: 80 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

MUM1, also named as EXPAND1, is an architectural component of the chromatin, which in response to DNA damage serves as an accessory factor to promote cell survival. It is a 53BP1-associated DNA damage-responsive factor. MUM1/EXPAND1 is different to another protein IRF4/MUM1.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87147-1-RR (MUM1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87147-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87147-1-RR against Human MUM1 were performed. The affinity constant is 0.974 nM.