

For Research Use Only

Mouse VPS13C Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87162-3-PBS



Basic Information

Catalog Number: 87162-3-PBS	GenBank Accession Number: NM_177184	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 320528	CloneNo.: 252231B7
Isotype: IgG	UNIPROT ID: Q8BX70	
Immunogen Catalog Number: AG29434	Full Name: vacuolar protein sorting 13C (yeast)	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean VPS13C concentration was determined to be 6.19 ng/mL in 3T3-L1 cell extract based on a 1.3 mg/mL extract load.



Sandwich ELISA standard curve of MP02898-1, Mouse VPS13C Recombinant Matched Antibody Pair - PBS only. 87162-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag29434. 87162-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL