

For Research Use Only

BLNK Recombinant monoclonal antibody, PBS Only

Catalog Number: 87170-3-PBS



Basic Information

Catalog Number: 87170-3-PBS	GenBank Accession Number: NM_013314.4	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 29760	CloneNo.: 252383A9
Isotype: IgG	UNIPROT ID: Q8WV28-1	
Immunogen Catalog Number: EG6610	Full Name: B-cell linker	
	Calculated MW: 50 kDa	
	Observed MW: 68-70 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

BLNK is a cytoplasmic linker or adaptor protein that plays a critical role in B cell development. This protein bridges B cell receptor-associated kinase activation with downstream signaling pathways, thereby affecting various biological functions. The phosphorylation of five tyrosine residues is necessary for this protein to nucleate distinct signaling effectors following B cell receptor activation. Mutations in this gene cause hypoglobulinemia and absent B cells, a disease in which the pro- to pre-B-cell transition is developmentally blocked. Deficiency in this protein has also been shown in some cases of pre-B acute lymphoblastic leukemia.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 87170-3-RR (BLNK antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 87170-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 87170-3-RR against Human LY86 were performed. The affinity constant is 1.13 nM.