

For Research Use Only

# BLNK Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87170-4-PBS



## Basic Information

Catalog Number:

87170-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG6610

GenBank Accession Number:

NM\_013314.4

GeneID (NCBI):

29760

UNIPROT ID:

Q8WV28-1

Full Name:

B-cell linker

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

252383A5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean BLNK concentration was determined to be 88.20 ng/mL in Ramos cell extract based on a 2.3 mg/mL extract load, 30.60 ng/mL in Raji cell extract based on a 1.2 mg/mL extract load and 6.50 ng/mL in A431 cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02893-2, Human BLNK Recombinant Matched Antibody Pair - PBS only. 87170-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6610. 87170-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL