

For Research Use Only

SLC9A3R1 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87174-2-PBS



Basic Information

Catalog Number: 87174-2-PBS	GenBank Accession Number: NM_004252	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 9368	CloneNo.: 252465F6
Isotype: IgG	UNIPROT ID: O14745	
Immunogen Catalog Number: AG30878	Full Name: solute carrier family 9 (sodium/hydrogen exchanger), member 3 regulator 1	
	Calculated MW: 39 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean EBP50/NHERF-1 concentration was determined to be 2219.3 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load, 162.0 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02896-1, EBP50/NHERF-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87174-2-PBS. Detection antibody: 87174-1-PBS. Standard: Ag30878. Range: 0.312-40 ng/mL.