

For Research Use Only

FAM84B Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:87176-1-PBS



Basic Information

Catalog Number: 87176-1-PBS	GenBank Accession Number: BC052957	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 157638	CloneNo.: 252473F4
Isotype: IgG	UNIPROT ID: Q96KN1	
Immunogen Catalog Number: AG13243	Full Name: family with sequence similarity 84, member B	
	Calculated MW: 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FAM84B concentration was determined to be 45.1 ng/mL in MDA-MB-453 cell extract based on a 1.8 mg/mL extract load, 4.3 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02903-1, FAM84B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87176-2-PBS. Detection antibody: 87176-1-PBS. Standard: Ag13243. Range: 0.156-20 ng/mL.