

For Research Use Only

UBD Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87177-1-PBS



Basic Information

Catalog Number:

87177-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3680

GenBank Accession Number:

BC012472

GeneID (NCBI):

10537

UNIPROT ID:

O15205

Full Name:

ubiquitin D

Calculated MW:

165 aa, 18 kDa

Purification Method:

Protein A purification

CloneNo.:

251913A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FAT10 concentration was determined to be 1700.5 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 805.0 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02908-1, FAT10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87177-2-PBS. Detection antibody: 87177-1-PBS. Standard: Ag3680. Range: 0.312-40 ng/mL.