

For Research Use Only

Serum amyloid A-1/SAA1 Recombinant monoclonal antibody, PBS Only (Detector)



Catalog Number: 87184-1-PBS

Basic Information

Catalog Number: 87184-1-PBS	GenBank Accession Number: GeneID (NCBI): 6288	Purification Method: Protein A purification
Source: Rabbit	UNIPROT ID: PODJ18	CloneNo.: 251324A8
Isotype: IgG	Full Name: serum amyloid A1	
Immunogen Catalog Number: EG4687	Calculated MW: 14kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Serum Amyloid A concentration was determined to be 5.6 ug/mL in human plasma.



Cytometric bead array standard curve of MP02910-1, Serum Amyloid A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87185-1-PBS. Detection antibody: 87184-1-PBS. Standard: Eg4687. Range: 0.625-80 ng/mL.