

For Research Use Only

PIGR Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 87186-2-PBS



Basic Information

Catalog Number: 87186-2-PBS	GenBank Accession Number: NM_002644.3	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5284	CloneNo.: 252416A1
Isotype: IgG	UNIPROT ID: P01833	
Immunogen Catalog Number: EG6497	Full Name: polymeric immunoglobulin receptor	
	Calculated MW: 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PIGR concentration was determined to be 71.1 ng/mL in human plasma.



Cytometric bead array standard curve of MP02912-1, PIGR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 87186-2-PBS. Detection antibody: 87186-1-PBS. Standard: Eg6497. Range: 0.156-20 ng/mL.