

For Research Use Only

VSIG1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 87189-1-PBS



Basic Information

Catalog Number: 87189-1-PBS	GenBank Accession Number: BC043216	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 340547	CloneNo.: 252228F5
Isotype: IgG	UNIPROT ID: Q86XK7	
Immunogen Catalog Number: EG6740	Full Name: V-set and immunoglobulin domain containing 1	
	Calculated MW: 387 aa, 42 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02900-1, Human VSIG1 Recombinant Matched Antibody Pair - PBS only. 87189-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg6740. 87189-1-PBS was HRP conjugated as the detection antibody. Range: 6.25-400 ng/mL.